



NIV masks and interchangeable elbow portfolio

## BiTrac NIV masks

The BiTrac NIV full face and nasal masks are single patient use masks that provide compatible ventilation for spontaneously breathing patients requiring noninvasive ventilation (NIV).

### Benefits for you and your patients



The mask provides a double-lip cushion.



The OmniClip is an adjustable piece on the forehead pad.



The headgear has adjustable straps with clips.



# BiTrac Select NIV masks interchangeable elbow portfolio

The new interchangeable elbow concept for full face and MaxShield masks



| Interchangeable elbow  | PN       | QTY |
|--|----------|-----|
| Standard - Used with a coaxial dual limb circuit on a Hamilton Medical ventilator.   | 10161785 | 10  |
| Anti-Asphyxia with leak* - Used for BiPAP devices without an exhalation port in the circuit. Compatible with Aerogen nebulizers.                           | 10161779 | 10  |
| Anti-Asphyxia* - Used for BiPAP devices with an exhalation port in the circuit.  | 10161784 | 10  |
| Anti-Asphyxia Nebulizer* - Used to administer nebulized medication on a BiPAP device. Exhalation port in the circuit. Compatible with Aerogen® nebulizers. | 10161780 | 10  |
| Standard Nebulizer - Used to administer nebulized medication on a Hamilton ventilator with a dual limb circuit. Compatible with Aerogen® nebulizers.       | 10161781 | 10  |
| Anti-Asphyxia Bronchoscopy* - Used for bronchoscopy on a BiPAP device.   | 10161782 | 10  |
| Standard Bronchoscopy - Used for bronchoscopy on a Hamilton ventilator.  | 10161783 | 10  |

\*Not to be used with ventilators

# BiTrac NIV masks

With standard interchangeable elbow for all Hamilton Medical ventilators, single patient use



## BiTrac® Select MaxShield

- Interchangeable standard elbow 22-mm ID
- Bonnet style headgear
- Recommended for patients 1 year or older and > 7 kg

| Size | Part number | QTY |
|------|-------------|-----|
| XXS  | 10161830    | 10  |
| XS   | 10161831    | 10  |

\*Only available in the US market.



## BiTrac NIV nasal mask\*

- Standard fixed elbow 22-mm ID
- OmniClip with silicone forehead pad and head strap
- Adult patients

| Size | Part number | QTY |
|------|-------------|-----|
| S    | 10161579    | 10  |
| M    | 10161578    | 10  |
| L    | 10161577    | 10  |
| XL   | 10161576    | 10  |



## BiTrac® Select NIV full face mask

- Interchangeable standard elbow 22-mm ID
- OmniClip with silicone forehead pad and head strap
- Adult patients

| Size | Part number | QTY |
|------|-------------|-----|
| S    | 10161769    | 10  |
| M    | 10161770    | 10  |
| L    | 10161771    | 10  |
| XL   | 10161772    | 10  |



## BiTrac® Select MaxShield

- Interchangeable standard elbow 22-mm ID and head strap
- Adult patients

| Size | Part number | QTY |
|------|-------------|-----|
| S    | 10161778    | 10  |
| L    | 10161777    | 10  |
| XL   | 10161776    | 10  |

# BiTrac NIV masks and circuit

With anti-asphyxia vented elbow for BiLEVEL and CPAP application, single patient use



## BiTrac® Select MaxShield\*

- Anti-asphyxia elbow 22-mm OD
- Bonnet style headgear
- Recommended for patients 1 year or older and > 7 kg

| Size | Part number | QTY |
|------|-------------|-----|
| XXS  | 10161832    | 10  |
| XS   | 10161833    | 10  |



## BiTrac NIV full face mask\*

- Standard fixed anti-asphyxia elbow 22-mm ID
- OmniClip with silicone forehead pad and head strap
- Adult patients

| Size | Part number | QTY |
|------|-------------|-----|
| S    | 10161765    | 10  |
| M    | 10161766    | 10  |
| L    | 10161767    | 10  |
| XL   | 10161768    | 10  |



## BiTrac® Select MaxShield\*

- Interchangeable anti-asphyxia elbow 22-mm ID and head strap
- Adult patients

| Size | Part number | QTY |
|------|-------------|-----|
| S    | 10161775    | 10  |
| L    | 10161774    | 10  |
| XL   | 10161773    | 10  |



## BiTrac NIV circuit

- 72" smooth ID tubing, swivel connector with exhalation port
- Available with or without attached filter

| Part number             | QTY |
|-------------------------|-----|
| 313-9075 with filter    | 10  |
| 313-9074 without filter | 10  |

\*Not for use with Hamilton Medical ventilators.



For more information and instructions for use, visit our e-Catalog at:  
[www.hamilton-medical.com/e-catalog](http://www.hamilton-medical.com/e-catalog)



**HAMILTON**  
**MEDICAL**

Distributed by

Hamilton Medical, Inc.

201 Edison Way, Unit A, Reno, NV 89502-2305

+1 800 426 6331

[info@hamilton-medical.com](mailto:info@hamilton-medical.com)

[www.hamilton-medical.com](http://www.hamilton-medical.com)

**Pulmodyne**<sup>®</sup>

Manufactured by

Pulmodyne, Inc

2055 Executive Dr, 46241 Indianapolis, US

+1 317 900 4093

[info@pulmodyne.com](mailto:info@pulmodyne.com)

[www.pulmodyne.com](http://www.pulmodyne.com)

ELO20210905N.01

The products shown here are not available for purchase by the general public. The information provided here is intended for healthcare professionals only. Always read the labels and follow the product's instructions for use. Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. All images are used for illustrative purposes only and may not accurately represent the product or its use. For all proprietary trademarks (®) and third-party trademarks used by Hamilton Medical AG see [www.hamilton-medical.com/trademarks](http://www.hamilton-medical.com/trademarks). ©2024 Hamilton Medical AG. All rights reserved.

NIV masks